


CZ BioHub_Take Home Test



Bo-Yen Chang

For Associated R&D Engineer

For Each Exercise:

- **Work**
- **Choice Made**


- **Challenges**
- **Solution**

Exercise 1

Exercise 1 — Explain the Work, Choice Made

Acquiring a multi-dimensional dataset using microscopy simulation

Work Explanation



```
z_step = 0.1
num_z_slices = 12
spacing = 500
channels = ['DAPI', 'FITC']
```

Choice made



- **pycromanager** — complex multi-dimensional acquisition and real-time image processing
- **pymmcore-plus** — lightweight, direct hardware control, without UI or advanced features.

Multi_dimensional_Acq_Events

```
# Step 6: Generate acquisition events, Using multi_d_acquisition_events
events = multi_d_acquisition_events(
    num_time_points=1,      # Set as 1, acquisit once at each event
    z_start=0, z_end=z_step*(num_z_slices-1.5), z_step=z_step, #Bas
    channels=channels,
    xy_positions=positions,
    order='tpcz'            # Acquisition order: time, position, chan
)
```

Exercise 1 — Challenges and solutions

Challenges

- Connect With the Micro-Manager
- Data Acquisition on Micro-Manager
- Dealing multi-dimensional data

Solutions

→ Look into the Official Document

- Connect to port 4827
- `pip install pycromanager`
- `import pycromanager`

→ Official Document

```
with Acquisition(directory='/path/to/saving/dir', name='acquisition_name') as acq:
    events = multi_d_acquisition_events(
        num_time_points=4, time_interval_s=0,
        channel_group='Channel', channels=['DAPI', 'FITC'],
        z_start=0, z_end=6, z_step=0.4,
        order='tcz')
    acq.acquire(events)
```

→ Two Method:

1. `events = multi_d_acquisition_events()`
2. Loop in the Acquisition Function

Exercise 1 — Verification Code to Check the Shape

```
def verify_z_step(acquisition_script_path, expected_z_step):
    """
    Verify that the z-step defined in the acquisition script matches the expected z-step.

    :param acquisition_script_path: The path to the acquisition script.
    :param expected_z_step: The expected z-step in micrometers.
    :return: None
    """
    # This is a simplified example and assumes you have a way to extract the z_step from the script.
    # In reality, you may need to parse the script or otherwise determine the z_step.
    z_step_defined_in_script = 0.1 # Example value

    if np.isclose(z_step_defined_in_script, expected_z_step):
        print("Z-step check passed.")
    else:
        print(f"Z-step mismatch: Expected {expected_z_step} um, got {z_step_defined_in_script} um")

# Load the dataset
path_to_dataset = r'C:\Users\chang\Desktop\2024_AssocRDEng_TakeHome\Exercise1\my_acquisition_6\my_acquisition_NDTiffStack.tif'
images = tiff.imread(path_to_dataset)
print(images.shape)

# Example usage:
expected_dims = (9, 12, 512, 512) # Adjust based on actual expected dimensions
expected_channels = ['DAPI', 'FITC']
expected_z_step = 0.1 # in micrometers

# Verify dimensions
if images.shape[:] != expected_dims[:]:
    print(f"Dimension mismatch: Expected {expected_dims[:]}, got {images.shape[:]}")
else:
    print(images.shape[1:])
    print("Dimension check passed.")
```

```
events = multi_d_acquisition_events(
    z_start=0,
    z_end=z_step*(num_z_slices-1.5),
    z_step=z_step,
)
```


```
# Based on the verification adjustment,
# '-1.5' could give in 12 slices
```

Exercise 2

Exercise 2 — Explain the Work, Choice Made


Modifying images in real-time during acquisition based on pixel value

Work Explanation



Choice made

Skimage.draw.disk



Process Image After Acquisition → `image_process_fn`

```
def img_process_fn(image, metadata):  
  
    ### send image and metadata somewhere ###  
  
    # this acquisition won't show a viewer or save data  
    with Acquisition(image_process_fn=img_process_fn) as acq:  
        ### acquire some stuff ###
```


Exercise 2 — Challenges and solutions

Challenges

- Cannot do disk drawing in Acquisition Function()
 - Acquisition Func() can only do event acquiring
- Data Type — `event`, `image` wrong
 - `image` should be 'numpy.array' but show 'dict'
 - `event` should be dict but show 'numpy.array' or 'queue'


Solutions

→ Have to use `image_process_fn` to call out def()

- Define two function
 - `process_image(image)`
 - Draw disk when image =700
 - `image_process_fn(image, metadata)`
 - Processing modified image (Save pics)
 - `metadata` is the event in this function

→ 1. Print in each process to see its 'type', 'shape'
2. Look into official document for function usage
3. Ask for micro-manager community

- `image_process_fn(image, metadata)`
 - has strict order for sending back value
- The `event` in `image_process_fn()` called `metadata`




Exercise 3


Exercise 3 — Explain the Work, Choice Made

Segmenting cell nuclei and measuring their shape deviations


Work Explanation



Segmentation




Eccentricities



Choice made

Segmentation - Cellpose:


- easy to access and specific to cell
- model()



Eccentricities - Skimage:



`measure.regionprops(mask).eccentricity`



Exercise 3 — Challenges and solutions

Challenges

- The images cannot be processed directly
- The data including 17 images, all of the images have multiple dimension

- `masks = model.eval()`

'list' object has no attribute 'ndim'

```
python Copy code  
  
# Process the images with Cellpose  
masks = model.eval(images, diameter=None, flow_threshold=None, channels=[0,0])
```

Solutions

→ To ensure 'image' was type: 'numpy array'

- using `io` in `scikit-image` to read the TIFF.

→ 1. Using 'for' loop to go over each images
2. Ensure images is the list of 2D arrays

- Only take first channel for processing

```
# Ensure images is a list of 2D arrays  
if images.ndim == 3:  
    ## If the image is 3D (multiple channels or Z-stacks), select  
    images = [img[0] if img.ndim == 3 else img for img in images]
```

→ `model.eval()`

- `model.eval()` must send back 4 variables
- `masks, flows, styles, diams = model.eval()`

What I have learned

- Learning the automation of microscopy from scratch
- Ask 'Official Document' and 'Community' for help
- Seeking for multiple tools to try more in different angles

Thanks!